

# FILT'RAY 3-month in-line filter

Ref. AT030743

Non-sterile



## DESCRIPTION

**FILT'RAY 3-month in-line filter**

**Ref. AT030743**

AQUATOOLS anti-Legionella and all-germ in-line filter.  
Suitable for filtering water to supply equipment used in the treatment of medical devices (e.g. endoscopes).  
Device made from MBS.  
Membrane type: tubular microfiltration.  
Membrane made from polyethylene.  
Pore cut-off threshold: 0.1µm nominal and 0.2µm absolute sterilising grade.  
Filtration flow rate: 6.3 lpm at 3 bar and 2.2 lpm at 1 bar (filter only, excludes any flow rate restriction in the mixer/tap).  
Maximum upstream pressure at the point-of-use: 5 bar.  
Maximum lifespan: 3 months from installation.  
Withstands thermal shocks: 70°C for 30 minutes.  
Withstands chemical shocks: active chlorine levels up to 100 ppm for 3 hours (cumulative total).  
Non-sterile filter.  
FILT'RAY 2G non-sterile disposable water filters are:  
CE marked Class I Medical Devices in Europe that comply with directive 93/42/CEE;  
Class II Medical Devices in the U.S.A. with FDA-510(K) clearance;  
Class I Medical Devices in Canada that comply with SOR/98-282.  
Packaged in boxes of 10, each filter is individually wrapped.

## ADVANTAGES



Anti-Legionella and all-germ



Lifespan: 3 months from installation



Bubl'air Wash: longer filter life



Rapid installation with no tools required

## TECHNICAL CHARACTERISTICS

**FILT'RAY 3-month in-line filter - Ref. AT030743**

Technology	Point-of-use filtration
Flow rate	4.5 lpm at 3 bar and 2.2 lpm at 1 bar
Norms	  

